

MANI DISPOSABLE SKIN STAPLER

MR MR conditional STERILE EO Sterilized using ethylene oxide 6pcs./Box


Product Name	Item	Order No.	Description
Manipler	Manipler S-2 35W	S235W	Wide size staple 35 count, single use
Manipler AZ	Manipler AZ 35W	AZ35WS	
	Manipler DSX 35W	DS35WS	

Non-sterile Non-Available in EU 1pc./Box

Product Name	Order No.	Description
Manipler Remover	SR-1S	Reusable skin staple remover


How to use

Manipler



Manipler AZ

The illustration shows the DSX type, but the usage is the same for the AZ type.



CE 0197

MANI, INC.

8-3 KIYOHARA INDUSTRIAL PARK, UTSUNOMIYA,
TOCHIGI, 321-3231, JAPAN
Phone: +81-28-667-7565 / E-mail: seg-ovs@mani.inc /
URL: www.mani.co.jp



Note : Design and specifications are subject to change without notice.

S0025



The best quality in the world,
to the world.


DISPOSABLE
SKIN STAPLER



Manipler AZ

Manipler AZ 35W

- User-friendly design**
Light grip force and light weight reduce user's stress and enable precise stapling
- Accuracy**
High strength steel material and practical staple shape ensure accurate wound closure
- Consistency**
High quality parts and structure achieve consistent performance



Manipler DSX 35W


- Light grip force**
Lighter grip force contributes to precise stapling and increases usability.
- Compact design**
Provides wider visibility.
Reduce medical waste at hospitals.
- 6-steps ratchet function**
Prevents jamming that is caused by lever return.
Increases easiness in staple positioning.



Manipler

Manipler S-2 35W

- Precision**
Ergonomic shape enables precise stapling as easy as possible
- Consistency**
High quality parts and structure achieve consistent performance
- Economical**
Compact design reduces warehouse space and disposal cost of medical waste



Manipler Remover

- Durability**
Durable structure ensures reliable performance



Not Available in EU

